## **Specifications TV4000LT**

## Tamson Visibility Bath 40 Litres Low Temperature - ASTM D445



Top lid stainless steel Integrated cooling High precision better than ± 0.02°C Large windows, no condensation Easy to operate RS232 communication

### General

The TV4000LT is a unique visibility bath with wide dimensions.

An integrated cooling system is able to lower the bath temperature down to minus 40°C within approximately two hours. The cooling power can be controlled over four stages so that the heat removal capacity is trimmed to the desired set point temperature, saving up to 50% of energy.

Due to its wide temperature range the bath can be used for multiple purposes varying from sub-zero temperature viscosity determination to calibration of sensors or thermometers. The large windows show the bath contents clearly. The windows are heated preventing build-up of condensate.

The bath can be emptied via a drain. If bath fluid expands at higher set point temperatures overflow of bath is protected via a drain. A LED fluid level indicator alarms when fluid level is too low.

Item	Unit	TV4000LT	
P/N 230V/50Hz		00T0460	
P/N 230V/60Hz		On request	
P/N 115V/60Hz		On request	
Power	[kW]	3 Max	
Used materials inside bath		Stainless steel, brass bearings	
Range		-40 +100°C -40 +212°F	
Reading		Standard °C, on request °F	
Setting	[°]	0.01	
Stability ±	[°C]	Better than 0.02	
Heating	[kW]	1.5 + 0.5 (2 heaters)	
Bath volume	[L]	40	
Opening	[mm]	260 x 240	
Depth bath	[mm]	300	
Length	[mm]	560	
Width	[mm]	492	
Height	[mm]	1150	
Noise level	[db]	61	
Weight	[Kg]	100	
CE	Conforms to CE regulation		

### Span

The bath temperature can be set from  $-40^{\circ}$ C to a maximum of  $+100^{\circ}$ C in steps of  $0.01^{\circ}$ C.

### Accuracy

Overall accuracy is better than  $\pm 0.02^{\circ}$ C over the whole temperature range. Ambient temperature drift <0.003°C/°C. The homogeneity is better than 0.02°C. Drift, accuracy and homogeneousness were measured with methanol as bath medium and are true min-max values found over eight hours.

#### Safety


The bath conforms to CE regulation. Further it is equipped with multiple safety features in order to guarantee safe around the clock operation. A mechanical safety thermostat protects the bath against high temperatures and will automatically reset itself after switching the bath off and on again. In case of over-temperature electronics and motor are switchedoff. The pump motor is protected against mechanical overload via a resettable fuse.




## **Specifications TV4000LT**

Tamson Visibility Bath 40 Litres Low Temperature - ASTM D445

## Cooldown of TV4000LT (methanol used as bath liquid)



### Temperature stability (methanol used as bath liquid)



PMT) Contact: G-Labo Germany 🕿 +49 6209 797100 😰 info@g-labo.de 🏠 www.g-labo.de

## Specifications TV4000LT

Standard cover and optional covers

TV4000LT is standard included with:				
P/N	Picture	Description		
0070400	2888	Cover with 7 openings: - 7 x ø51 mm opening - 2 x ø12.5mm opening for thermometer		
23T2400	$\bigcirc$	7 * lid for ø 51 mm opening		

	Optional covers for TV4000LT:					
	P/N	Picture	Description			
	23T2401 -	00000	Cover with 8 openings: - 8 x ø51 mm opening - 2 x ø12.5mm opening for thermometer			
		$\bigcirc$	8 * lid for ø 51 mm opening			
	23T2402 -		Cover with 8 openings: - 8 x ø60 mm opening - 2 x ø12.5mm opening for thermometer			
		$\bigcirc$	8 * lid for ø 60 mm opening			
	23T2403		Cover with 7 openings: - 4 x ø51 mm opening - 3 x ø60 mm opening - 2 x ø12.5mm opening for thermometer			
		$\bigcirc$	4 * lid for ø 51 mm opening			
		$\bigcirc$	3 * lid for ø 60 mm opening			

## Specifications TV4000LT

Accessories

Accessories						
P/N	Picture	Description				
19T4024		E20 Digital contact thermometer.Two decimal readout. Temperature range -40°C to 140°C. Please see specification sheet "E20 Thermometer" for more information.				
14T0303		Adapter to insert an E20 thermometer in the opening of the cover				
Viscosity Accessories		Please see specification sheet "Viscosity Accessories"				