# Specifications TC16 - TC40 - TC58

## Thermostatic Bath and Circulator



#### General

The TC16, TC40, and TC58 are used for precise temperature control. The baths offer a choice in different bath volumes and offer high reliability for use under rough circumstances.

#### Accuracy

The insulation of the bath and electronic design result in a stable working temperature of  $\pm 0.02^{\circ}$ C for the TC16, TC40, and TC58. The stability was measured in water at 50°C. The set point temperature can be set in steps of 0.1°C in the range of 0°C up to 250°C (32°F..482°F). The accuracy on the display is displayed in 0.1°C. However the controller has an internal accuracy of 0.01°C. This two decimal-readout also is available using the RS232 communication port and the free Tamcom software.

### **Temperature readout**

Standard available in °C, on request in °F.

## Pump

Three pumps can be offered:

Suction/pressure pumpPerformance pump

- Standard pump
- 00T0253 24T0399

Item	Unit	TC16	TC40	TC58		
P/N 230V/50-60Hz		00T0671	00T0681	00T0691		
P/N 115V/60HZ		00T0861	00T0851	00T0881		
Power	[Watt]	1600	3000	3000		
Materials	Used inside bath: stainless steel 304, brass					
Range	Ambient 250°C/482°F					
Reading	Standard °C, °F on request					
Setting	[°]	0.1				
Stability ±	[°C]	0.02	0.02	0.02		
Heating	[W]	1400**	2800	2800		
Bath volume	[L]	16	40	58		
Pump	[mBar]	Max pressure 300*				
Pump	[L/min]	Max flow 7.5*				
Opening	[mm]	180*210	420*275	420*275		
Depth bath	[mm]	220	200	300		
Length	[mm]	480	705	705		
Width	[mm]	295	375	375		
Height	[mm]	480	440	590		
Weight	[Kg]	19	26.8	35		
CE	Conforms to CE regulation					
* Optional pump with pressure max 1Bar, flow max 16L/min,						

\* Optional pump with pressure max 1Bar, flow max 16L/min, reference: 24T0399. See specification sheet "pump graph" \*\* Additional heater 1400W optional to heat larger loads when circulating

#### Span

The baths can be operated from ambient up to  $250^{\circ}$ C (482°F). Due to the "self-heating" of the stirring mechanism (friction-heat) and when using a top lid, the minimum temperature lies approximately 5°C above the ambient temperature. When using the standard built-in cooling coil the minimum operating temperature lies approximately  $\pm$  5°C above the temperature of the cooling liquid circulating inside the cooling coil. Minimum working temperature when using external cooling lies around 10°C. Lower temperatures generate problems with built up of condensate in the inner insulation of the bath.

#### Safety

The bath conforms to CE regulation. It is further equipped with a mechanical resettable safety thermostat.

#### Optional

See table on the next page.

est 1878

MT Contact: G-Labo Germany 🏗 + 49 6209 797100 🛕 info@g-labo.de 🏠 www.g-labo.de

# Specifications TC16 - TC40 - TC58

## Thermostatic Bath and Circulator

#### Operating

- 1 LED Temperature display
- 2 Item
- 3 Sequence
- 4 Up & Down
- 5 Over temperature safety
- 6 Power

The thermostatic bath is easy to operate. The four keys will set all the parameters. Using the up and down arrow the desired working temperature or SP can be set in steps of 0.1°C. The bath temperature is displayed with a resolution of  $0.1^{\circ}$ C.



	Ite	m	TC16	TC40	TC58		
	Remote PT100		On request				
	RS232 serial communica	ation	02T3025				
	Software "Tamcom"		Free of charge				
	Suction/pressure pump		00T0253				
	Performance pump		24T0399				
	Bath drain		Standard included				
	Levelling platform		07T0210	07T0211	07T0212		
C	Lid levelling platform	A REAL PROPERTY AND A REAL	Included in 07T0210	03T2270			
2	Gabled lid			03T2210 (In combination with 07T0211/07T0212)			
	Top lid with cover of one set of nine rings		00T0254				
	Top lid with cover of two sets of six rings (opening without rings ø180 mm)			00T(	0255		
	Top lid with cover of six sets of three rings (opening without rings ø95 mm)			00T(	0256		



# Tamson Instruments Specification sheet

# Specifications TC16 - TC40 - TC58

# Thermostatic Bath and Circulator

Item		TC16	TC40	TC58	
Cooling circulator TLC10-3 230V/50Hz		00T0050			
Cooling circulator TLC10-3 230V/60Hz			00T0051		
Cooling circulator TLC10-3 115V/60Hz		00T0052			
Cooling circulator TLC15-5 230V/50Hz			<b>00</b> T	0565	
Cooling circulator TLC15-5 230V/60Hz			00T	0567	
Cooling circulator TLC15-5 115V/60Hz			00T	0570	
Additional heater (1400W)		<b>25T0194</b> - 230V			

**TC40**