

ASTM D 381 - ISO 6246 - IP 131 - DIN 51784

SCOPE

Standard Test Method for Gum Content in Fuels by Jet Evaporation.

EXISTENT GUM IN FUELS BY JET EVAPORATION

- **AIR/STEAM JET CLASSIC** -

- **REF 941303** -

MAIN CHARACTERISTICS

Air and steam jet evaporation, with 5 wells, including:

- ✓ Electronic temperature control with digital display and flowmeter +155°C +/-0,5°C et +232°C +/-2°C
- ✓ Built-in steam super heater
- ✓ Thermometer 3C
- ✓ Beakers
- ✓ Conical outlets (5 pieces)

SCOPE OF DELIVERY

AIR/STEAM JET CLASSIC delivered ready for use with:

REF 9413031 Super heater for 941303 apparatus

CONSUMABLES

REF 12875 Filter funnel, sintered glass, 125 ml capacity
REF 16135 Ungraduated beaker without Neck 50 x 78 MM, (150 ml)
REF 27113 50 ml borosilicate glass cylinder

OPTIONS ACCESSORIES

REF 9413080 Air compressor 2 stages suitable to deliver correct flow of air, 300 l capacity, for ref. 941302 and 941303. For use on AC 400 V - 3 ph - 7 A - (1670x600x760mm) ± 125 kg.

REF 941308 Steam generator (400 V - 3 phases) for use with ref. 941303 apparatus, equipped with sensitive constant level device to maintain vapor output at requested quantity, supplied with flexible connecting pipe, output capacity : 2 to 12 kg, pressure 0-100 psi, with automatic security control, loading 12 Kw
 For use on AC 400 V - 3 ph - 20 A - (750x530x920mm ±75 kg).

REF 941307 Overpressure pump for steam generator (230 V) for ref. 941308

REF 11491 ASTM thermometer (3C)

REF 9417444 Pair of gloves



**AIR/STEAM JET CLASSIC
REF. 941303**



REF. 941308

ORDERING INFORMATION

REF 941303
AC 230 V - 50 Hz - 16 A

(770x380x570mm) ±60 kg.

CONTACT : sales@normalab.com

NORMALAB FRANCE SAS
 ZA Caux Multipôles 1 - F-76 190 Valliquerville
 Tel. : +33 232 700 100
 Fax : +33 232 704 732

www.normalab.com


REF 941303 Leaflet 15

DISTRIBUTED BY

G-Labo Germany
 Bgm.-Horneffstr.26
 69509 Mörlenbach
 Tel.: + 49 6209 797100
 Fax: + 49 6209 797101
 Mail: info@g-labo.de
 Web: www.g-labo.de



REF. 9413080



REF. 941301 & REF. 941304

CALIBRATION EQUIPEMENTS

REF 941301 Air calibration kit
AC 230 V - 0.3A - 50 Hz (680x340x180mm) ±5 kg
REF 941304 Vapour calibration kit

CONTACT : sales@normalab.com

NORMALAB FRANCE SAS
ZA Caux Multipôles 1 - F-76 190 Valliquerville
Tel. : +33 232 700 100
Fax : +33 232 704 732

www.normalab.com

REF 941303 Leaflet 15

DISTRIBUTED BY

G-Labo Germany
Bgm.-Horneffstr.26
69509 Mörlenbach
Tel.: + 49 6209 797100
Fax: + 49 6209 797101
Mail: info@g-labo.de
Web: www.g-labo.de