

ASTM D 482 - D 874

SCOPE

Standard Test Method for Ash from Petroleum Products.

ASH FROM PETROLEUM PRODUCTS - REF 23963 -

MAIN CHARACTERISTICS

- ✓ Muffle furnace model L3/11
- ✓ 3 liters capacity
- ✓ Maximum temperature of 1100°C.

SCOPE OF DELIVERY

ASH FROM PETROLEUM PRODUCTS delivered ready for use with :

- REF 23964** Set of heating elements
REF 23965 Thermocouple K (160mm)

NECESSARY ACCESSORIES (for ASTM D482)

- REF 10532** Porcelain dish, 100 mL, flat bottom (pack of 6)
REF 10687 Stainless curved crucible tong (pack of 2)
REF 10702 Glass balls diam. 2mm (1kg)
REF 10770 Porcelain plate for below dessicator diam. 200mm
REF 9411300 Platinum crucible (100mL)

NECESSARY ACCESSORIES (for D874)

- REF 10508** Porcelain crucible, 40mL tall form (pack of 6)
REF 10509 Porcelain crucible, 80mL tall form (pack of 6)
REF 10510 Porcelain crucible, 120mL tall form (pack of 6)
REF 10515 Porcelain crucible, 80mL low form (pack of 6)
REF 10532 Porcelain dish, 100mL, flat bottom (pack of 6)
REF 10687 Stainless curved crucible tong (pack of 2)
REF 10702 Glass balls diam. 2mm (1kg)
REF 10770 Porcelain plate for below dessicator
REF 27203 Soda lime glass dessicator diam. 200mm
REF 15377 Silica crucible (50 mL tall form)
REF 15386 Silica crucible (75 mL tall form)
REF 15387 Silica crucible (130 mL tall form)
REF 15671 Platinum crucible (50mL)



ORDERING INFORMATION

REF 941281
AC 230 V - 50 Hz - 6 A (\pm 20 kg)

Internal size : 160 x 140 x 100 mm
External size: 380 x 370 x 420 mm

CONTACT : sales@normalab.com

NORMALAB FRANCE SAS
ZA Caux Multipôles 1 - F-76 190 Valliquerville
Tel. : +33 232 700 100
Fax : +33 232 704 732

www.normalab.com

REF 23963 Leaflet 14rev1

DISTRIBUTED BY

G-Labo Germany
Bgm.-Horneffstr.26
69509 Mörlenbach
Tel.: + 49 6209 797100
Fax: + 49 6209 797101
Mail: info@g-labo.de
Web: www.g-labo.de