

ASTM D 1322 - ISO 3014 - IP 57 - DIN 51406

SCOPE

Standard Test Method for Smoke Point of Kerosene and Aviation Turbine Fuel

**DETERMINATION OF THE SMOKE OF KEROSENE AND AVIATION TURBINE FUEL
- REF 941680 -**

MAIN CHARACTERISTICS

- ✓ Smoke point apparatus ready for use
- ✓ interchangeable oil container

SCOPE OF DELIVERY

Ready to use with:

REF 9416801 Candle (oil container)

CONSUMABLES

REF 9416802 Pack of 12 wicks



ORDERING INFORMATION

REF 941680

(W) 200x (D) 200x (H) 400 mm - (±4 kg)

CONTACT : sales@normalab.com

NORMALAB FRANCE SAS
ZA Caux Multipôles 1 - F-76 190 Valliquerville
Tel. : +33 232 700 100
Fax : +33 232 704 732

www.normalab.com

Ref 941680 Leaflet 14rev1

DISTRIBUTED BY

G-Labo Germany
Bgm.-Horneffstr.26
69509 Mörlenbach
Tel.: + 49 6209 797100
Fax: + 49 6209 797101
Mail: info@g-labo.de
Web: www.g-labo.de