

**ASTM D 97 - D 2500 - ISO 3015 & ISO 3016 - IP 15 & IP 219**

**SCOPE**

Standard Test Method for Pour Point of Petroleum Products

**CLOUD AND POUR POINT BATH**  
**- CPP BASIC-**  
**- REF 941584**

**MAIN CHARACTERISTICS**

- ✓ Temperature bath with 5 wells
- ✓ Stirrer
- ✓ One set of 6 thermometers 5C
- ✓ Cooled by dry ice or cryostat

**SCOPE OF DELIVERY**

5 felt discs sold in packs of 10 under **REF 9415901**  
5 felt rings sold in packs of 10 under **REF 9415902**

**REF 19439** Test tube  
**REF 11492** ASTM thermometer (5C)

**NECESSARY ACCESSORIES**

**REF 941597** Dry Ice Machine  
**REF 23205** Cryostat TLC 80 11 litres (-80°C / +20°C)  
**REF 21131** Insulated pipe for cryostat internal Ø 7mm  
**REF 23204** Cryostat TLC40/14 liters (-45°C...+20°C ; stab ± 0.05K)  
**REF 11493** ASTM thermometer (6C)  
**REF 11530** ASTM thermometer (61C)

**REFERENCE FLUIDS**

**REF 60002** Cloud point (Diesel) ; -5,3°C  
**REF 60003** Pour Point (Diesel) ; -31,7°C



**ORDERING INFORMATION**

**REF 941584**  
For use on AC 230 V, 50 Hz - 1 A

(W) 290x (D) 360 x (H) 250 mm (± 9 kg)

CONTACT : [sales@normalab.com](mailto:sales@normalab.com)

NORMALAB FRANCE SAS  
ZA Caux Multipôles 1 - F-76 190 Valliquerville  
Tel. : +33 232 700 100  
Fax : +33 232 704 732

[www.normalab.com](http://www.normalab.com)

CPP BASIC Leaflet 14rev2

**DISTRIBUTED BY**

G-Labo Germany  
Bgm.-Horneffstr.26  
69509 Mörlenbach  
Tel.: + 49 6209 797100  
Fax: + 49 6209 797101  
Mail: [info@g-labo.de](mailto:info@g-labo.de)  
Web: [www.g-labo.de](http://www.g-labo.de)