### MODEL CPP CLASSIC

# COLD AND POUR POINT TEST CABINET



#### STANDARDS

ASTM D 97, D 2500, IP 15, ISO 3015 & 3016 and related methods.

#### SCOPE

The pour point test method covers and is intended for use on any petroleum product. The cloud point test method covers only petroleum products and biodiesel fuels that are transparent in layers 40 mm in thickness, and with a cloud point below 49°C.



#### **SPECIFICATIONS**

**The CPP test cabinet** is mechanically refrigerated by CFC free Gas. There are 4 compartments (0, -18, -33, -51°C), each with 4 wells. It comes with useful glassware and thermometers. The CPP Classic has an internal cooling system.

#### **Key features**

- 4 electronic temperature regulators
- Refrigerated by CFC free gases
- 4 compartments (0, -18, -33, -51°C) each with 4 wells
- Temperature stability: +/- 1°C
- Adjustment accuracy for each bath: 1°C.
- · Bath volume: 3 liters

#### ORDERING INFORMATION

## NORMALAB production - Half-automated model of cloud and pour point tester

#### 941592

#### Scope of delivery:

CPP Classic is delivered ready to use with:

- 4 tube holders
- 16 glass test jars (P//N 19439)
- 4 thermometers ASTM 5C \* (P/N 11492)
- 4 thermometers ASTM 6C \* (P/N 11493)
- 16 cork stoppers (P/N 0941400154)
- \* Without certificate, calibration on request

#### Site requirements:

- Power supply: 230 V 50/60 Hz 1815 W
- Dimension: (W) 910 x (D) 630 x (H) 1250 mm
- Weight: 190 kg

Fluid: 12 L ethanol.



**CONTACT**: sales@normalab.com

Normalab FRANCE SAS

ZA Caux Multipôles 1 - 76190 Valliquerville

Tel.: +33 232.700.100 Fax: +33 232.704.732

#### **DISTRIBUTED BY**

G-Labo Germany Bgm.-Horneffstr.26 69509 Mörlenbach

Tel.: + 49 6209 797100 Fax: + 49 6209 797101 Mail: info@g-labo.de Web: www.g-labo.de