

ASTM D 113 - NFT 66006 & DIN 52013

SCOPE

Standard Test Method for Ductility of Bituminous Materials.

DUCTILITY OF BITUMINOUS MATERIALS
-DM 150 CLASSIC -
-REF 9418642-

MAIN CHARACTERISTICS

- ✓ Ductility machine 3 posts
- ✓ Length 150 cm, with a circulation bath and thermostating coil
- ✓ Variable speed from 1 to 99 mm / min in steps of 0.1 mm,
- ✓ Background consisting of two plates of opaline reflective effect
- ✓ Draincock
- ✓ Drive motor with worm screw,
- ✓ Mobile assembly built in bearing
- ✓ Digital measurement display

SCOPE OF DELIVERY

DM 150 CLASSIC delivered ready for use with :

- REF 941868** Mould without brass plate
REF 941869 Brass plate
REF 941865 Opal glass for ductility meter 150 cm
REF 11532 ASTM thermometer (63C)

ACCESSORIES

- REF 23203** Cryostat TLC 30 - 5 litres capacity (-30 to +60°C / ±.1)
REF 9418671 Transparent cover for ductility meter (1m50)



ORDERING INFORMATION

REF 9418642
DM 150 CLASSIC For use on AC 230 V - 50 Hz - 3 A

(W) 2200x (D) 330x (H) 310 mm (± 85 kg)

CONTACT : sales@normalab.com

NORMALAB FRANCE SAS
ZA Caux Multipôles 1 - F-76 190 Valliquerville
Tel. : +33 232 700 100
Fax : +33 232 704 732

www.normalab.com

REF 9418642 Leaflet 14rev1

DISTRIBUTED BY

G-Labo Germany
Bgm.-Horneffstr.26
69509 Mörlenbach
Tel.: + 49 6209 797100
Fax: + 49 6209 797101
Mail: info@g-labo.de
Web: www.g-labo.de