

**ASTM D 721 - ISO 2908 - IP 158 - DIN 51571 & D 3235**

**SCOPE**

Standard Test Method for Oil Content of Petroleum Waxes.

**OIL CONTENT OF PETROLEUM WAXES**  
**- OCW CLASSIC -**  
**- REF 941795 -**

**MAIN CHARACTERISTICS**

- ✓ Complete apparatus with cooling bath 3 filter sticks
- ✓ Air pressure regulator
- ✓ Evaporation assembly 4 weighing bottles and thermometer
- ✓ Diver with ground sintered
- ✓ Glass door
- ✓ Manifold air distribution on rack, with 4 bottles of 15 ml and weigh thermometer

**SCOPE OF DELIVERY**

OCW CLASSIC delivered ready for use with :

**REF 23240** Weighing bottle (15 ml)

**REF 517564** Cooling sample tube

**REF 517565** Filter stick B24 x3

**CONSUMABLE**

**REF 517565** Filter stick B24 x3



**ORDERING INFORMATION**

**REF 941795**

OCW CLASSIC For use on AC 230 V - 50 Hz - 0.3 A

(W) 600x (D) 300x(H) 450 mm (± 15 kg)

CONTACT : [sales@normalab.com](mailto:sales@normalab.com)

NORMALAB FRANCE SAS  
ZA Caux Multipôles 1 - F-76 190 Valliquerville  
Tel. : +33 232 700 100  
Fax : +33 232 704 732

[www.normalab.com](http://www.normalab.com)

REF 941795 Leaflet 14rev1

DISTRIBUTED BY

G-Labo Germany  
Bgm.-Horneffstr.26  
69509 Mörlenbach  
Tel.: + 49 6209 797100  
Fax: + 49 6209 797101  
Mail: [info@g-labo.de](mailto:info@g-labo.de)  
Web: [www.g-labo.de](http://www.g-labo.de)