

ISO 1430 - ASTM in Progress

DETERMINATION OF PARTICLE POLARITY OF BITUMEN EMULSIONS

- **NEB 110 MODEL** -
- **REF 941890** -

MAIN CHARACTERISTICS

NEB 110 allows the determination of the sign of the electric charge carried by particles in ionic emulsions.
A direct or rectified current is passed through a bitumen emulsion between 2 parallel electrodes.

SCOPE OF DELIVERY

REF 12449 Graduated lipped beaker – 250 ml
REF 941891 Complete electrode with cover



CONTACT : sales@normalab.com

NORMALAB FRANCE SAS
ZA Caux Multipôles 1 - F-76 190 Valliquerville
Tel. : +33 232 700 100
Fax : +33 232 704 732

www.normalab.com

REF 941890 Leaflet 14rev1

DISTRIBUTED BY

G-Labo Germany
Bgm.-Horneffstr.26
69509 Mörlenbach
Tel.: + 49 6209 797100
Fax: + 49 6209 797101
Mail: info@g-labo.de
Web: www.g-labo.de